


LOGIN GUIDE


This login guide provides the specific login information for your device. Please refer to the manual for all information on using myDiabby Healthcare*.




- 1 **Create a Dexcom Clarity account** from this site:
<https://uam2.dexcom.com//index.html#/register>



- 2 **Download the Dexcom Clarity Uploader** from **this page**:
<https://clarity.dexcom.eu/?locale=fr-FR#/home-user-upload>
Then **upload the monitor data to Dexcom Clarity** by following the instructions on the screen.





- 3 **Log in to the myDiabby account** from a computer **OR** the mobile application.
- 4 **Click on the "My Connected Objects" tab** and then on **"Connect my continuous glucose monitor"**.



5 Choose the Dexcom button, then click on the **"CONNECT MY HEALTH APPLICATION"** button




6 Enter your **Dexcom account login details and password**. Read the terms and conditions of use, then **enter your name and click "Allow"** to accept.



7 This message appears on the screen and confirms that the **Dexcom account is connected to myDiabby**. Click on **"Back"** to return to myDiabby Healthcare.



Continuous blood glucose data is now shared on myDiabby Healthcare !



PLEASE NOTE

Upload the Dexcom meter data to Dexcom Clarity regularly so that the healthcare team can see the results on myDiabby.

> Once Dexcom and myDiabby are connected, **Dexcom automatically sends data to myDiabby once an hour**. The data displayed is 3 hours behind the last available data.

> **If you are using a Dexcom G6 in combination with the t:slim X2 Tandem Insulin Pump, you only need to upload the pump to myDiabby for the Tandem Pump + G6 monitor data to appear on myDiabby Healthcare.**

If you have any questions about the myDiabby application:
support@mydiabby.com or by phone

France : **01 76 40 01 78**

Belgique : **02 320 11 96**



MDHC SAS
66 avenue des Champs Elysées
75008 Paris
FRANCE

*The manual is available for download and consultation on the myDiabby Healthcare platform.



0459